

In-gel Digestion of Excised Protein Spots

All steps at room temperature unless noted otherwise. Adjust volumes as required to cover gel pieces for large bands/spots.

1. **Washing and Destaining Coomassie blue stained gels.**
 - a. Rinse gels 2 times for 15 min in ddH₂O.
 - b. Excise bands or spots of interest, cut into cubes (< 1mm) and place in low retention microfuge tubes (Siliconized, e.g. Fisher 02-681-311).
 - c. Cover pieces with 40 μ l 100 mM ammonium bicarbonate (ambic) plus 40 μ l acetonitrile and incubate 15 min.
 - d. Remove supernatant and cover again with 1:1 acetonitrile/ 100 mM ambic. Use a gel-loading tip to prevent sucking up gel pieces.
 - e. Repeat steps 1c and 1d until gel pieces are clear.
2. **Reduction and alkylation (Start here for silver stained gels).**
 - a. Add 40 μ l of 10 mM DTT in 100 mM ambic and incubate at 56°C for 45 min.
 - b. (optional) Remove supernatant, immediately add 40 μ l of 55 mM iodoacetamide in 100 mM ambic and incubate 30 min in the dark.
 - c. Remove supernatant and cover pieces with 40 μ l 100 mM ambic and incubate 5 min
 - d. Add 40 μ l of acetonitrile.
 - e. Remove supernatant and dry gel pieces completely in a speed vac.
3. **Digestion and extraction.**
 - a. Add 40 μ l of trypsin solution and incubate 45 min on ice. Add more solution if all the liquid is absorbed.
 - b. Remove excess solution, cover pieces with digestion buffer and incubate 16 h at 37°C
 - c. Remove supernatant and save in a low retention microfuge tube.
 - d. Add 20 μ l of 25 mM ambic and incubate 15 min.
 - e. Add 20 μ l of acetonitrile to make 50% solution and incubate 15 min.
 - f. Remove supernatant and add to solution from step 3c.
 - g. Add 20 μ l of 5% formic acid or 1% TFA to gel pieces and incubate 15 min.
 - h. Add 20 μ l of acetonitrile to make 50% solution and incubate 15 min
 - i. Remove supernatant and add to solution from step 3c.
 - j. Repeat steps 3g-3i
 - k. Add 10 mM DTT to give a final concentration of 1mM.
 - l. Dry extracted peptides in a speed vac.
 - m. Resuspend in 20 μ l of 1% TFA
4. **Zip-Tip Purification**
 - a. Wet Zip-Tip (Millipore ZTC18S096) with 10 μ l Wetting/Elution solution, expel. (wait for tip to fill and expel completely to avoid bubbles in the matrix)
 - b. Rinse 2X with wash solution. Expel and hold plunger down after 2nd wash.
 - c. Bind peptides to tip by filling and expelling resuspended peptides through zip-tip several times (more cycles for more complex mixtures). At the last cycle leave the solution in the tip to allow for preparation of a tube for elution.
 - d. Wash tip 2X with wash buffer. Expel and hold plunger after 2nd wash.
 - e. Elute peptides with 2 μ l of elution solution into the bottom of the fresh tube, flush the solution through the tip 2-3X and expel after the final rinse. Repeat elution step by suspending the droplet on the side of the tube. Combine with the first eluate. The sample is ready for MALDI-TOF analysis. The samples could be concentrated slightly by allowing the solvent to evaporate but this is usually not necessary.

Stock Solutions

- 1 M ammonium bicarbonate
0.79 g ammonium bicarbonate in 10 ml ddH₂O (pH should be ~7.8)
- 1 M DTT
0.15 g DTT in 1 ml ddH₂O
- 1 M CaCl₂
0.11 g CaCl₂ in 1 ml ddH₂O
- Trypsin, Sequencing grade
25 mg dissolved in 0.25 ml ddH₂O

Digestion Solutions

Prepare all solutions fresh before performing the digest.

100 mM Ammonium bicarbonate 100 µl 1M ammonium bicarbonate 900 µl ddH ₂ O	55 mM iodoacetamide in 100 mM ambic 10 mg iodoacetamide 100 µl 1M ammonium bicarbonate 900 µl ddH ₂ O
10 mM DTT in 100 mM ambic 10 µl 1 M DTT 100 µl 1M ammonium bicarbonate 900 µl ddH ₂ O	Trypsin solution (prepare enough to process the number of samples in each expt., 40 µl per sample) 34 µl/sample Digestion buffer 6 µl/sample 0.1 mg/ml Trypsin
Digestion buffer 50 µl 1M ambic 5 µl 1 M CaCl ₂ 945 µl ddH ₂ O	10 mM DTT 10 µl 1 M DTT 990 µl ddH ₂ O
25 mM ambic 25 µl 1 M ammonium bicarbonate 975 µl ddH ₂ O	1% TFA 10 µl Trifluoroacetic acid 990 µl ddH ₂ O
Zip-Tip wash solution 100 µl 1% TFA 900 µl ddH ₂ O	Zip-Tip wetting/elution solution 500 µl acetonitrile (HPLC or UV grade) 100 µl 1% TFA 400 µl ddH ₂ O